
Coapt

ENDOTINE Triple™ Device

Instructions for Use

English	1
Dansk	2
Deutsch	3
Español	5
Français	6
Italiano	8
Nederlands	10
Português	11
Svenska	13
Polski	14
Magyar	16
Ελληνικά	17
Русский	19


1


2


4

3

5

ENDOTINE Triple™ Device Instructions for Use

Available Models

Part Number	Tine Length	Part Name
CFD-010-0167	3.0 mm	ENDOTINE Triple™

DESCRIPTION

The ENDOTINE Triple™ device is a bioabsorbable fixation implant. This device is SINGLE USE ONLY and is supplied sterile. The device and its insertion tools are for SINGLE USE ONLY and are supplied sterile.

INDICATIONS

The ENDOTINE Triple™ is intended for use in subperiosteal browplasty surgery. The ENDOTINE Triple™ device is specifically indicated for use to fixate the subcutaneous tissues to the cranial bone in browplasty.

CONTRAINDICATIONS

1. Patients appearing to have very thin bones which might imply an inadequate cranial thickness for the ENDOTINE Triple™ device; the correctly drilled hole extends to a depth of $3.00\text{mm} \pm 0.025$.
2. Situations where internal fixation is otherwise contraindicated (e.g. infection).
3. Thin, atrophic scalp.
4. Any known allergy or foreign-body sensitivities to plastic biomaterial.

LABELING, PACKAGING and STERILIZATION

The ENDOTINE Triple™ device and insertion tools are sterilized by gamma irradiation. **Do NOT re-sterilize. The device should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of punctures, tampering, water contamination or other damage.**

STORAGE INSTRUCTIONS

1. Store at room temperature (15 to 24°C or 60 to 75°F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the label.

INSTRUCTIONS FOR USE

(For endoscopic or open techniques)

1. Prior to implant usage, verify intraoperatively that the scalp thickness is appropriate for the ENDOTINE Triple™ devices.
2. Perform the dissection in the usual fashion, making sure to achieve adequate release.
3. To minimize the risk of creating a full thickness cranial hole with the ENDOTINE Triple™ drill bit, verify that the hole position is medial to the temporal fusion line and anterior to the coronal suture (**Figure 1**). Avoid the sagittal midline.

4. If possible, plan the cranial hole placement such that the implant will ultimately lie under hair bearing tissue. Typically, for sagittal, paramedian incisions, the implant lies somewhat lateral and/or anterior to the incision. For coronal or transverse incisions, the implant lies anterior to the incision.
5. Prepare the cranial bone hole using the enclosed drill bit (Figure 2) applied perpendicular to the bone. **After drilling, check that no residual loose bone material is present in the drill hole.**
6. Remove the ENDOTINE Triple™ device and insertion tool from the package.
7. Carefully place the device through the scalp incisions to avoid dislodgement from the insertion tool. Then place the device post into the pre-drilled cranial hole with moderate pressure until the device platform is resting on the cranium. This should be verified by visual inspection (**Figure 3**).
8. Release the insertion tool from the implant by rotating the insertion tool in a counterclockwise direction (**Figure 4**). The platform of the ENDOTINE Triple™ device should remain flush against the cranium.
9. Perform the brow elevation and apply manual pressure to the scalp to ensure penetration of the tissue by the device (**Figure 5**).
10. If a drain is used, manipulate the subcutaneous portion such that it is maximal distance from the implant.
11. Use caution during skin closure to avoid elevating the scalp tissue off the implant.
12. A gentle pressure dressing may be appropriate upon completion of the case.

PRECAUTIONS

1. If possible, the implant should lie under hair bearing scalp at the conclusion of the procedure. Placement under forehead or alopecic skin may lead to unacceptable visibility of the device prior to absorption.
2. The patient is to be informed of possible adverse effects including palpability of the device prior to resorption, discomfort, and surgical risks. Tenderness and visible fullness over the device for several months has been noted, especially in patients with thin scalps. These symptoms and/or wound problems may lead to a decision to surgically remove the device prior to its full resorption.
3. Local cyst and/or sinus formation may occur and may be treated by aspiration or device removal.
4. High speed drilling (approximately 1000 rpm and greater) and repeat drilling may create holes that are wider than intended for a tight grip of the ENDOTINE Triple™ device. Single pass, low speed drilling reduces the occurrence of this issue.
5. Avoid any movement of the scalp, which may dislodge the tissue from the tines.
6. Support the head when turning or moving the patient.

WARNINGS

1. A drill bit is supplied for use in full size adult patients to create a hole in the cranium above (medial to) the temporal fusion line and anterior to the coronal suture. However, the surgeon should note any penetration beyond the inner table of the cranium and, if observed, should seek follow-up and neurosurgical consultation as appropriate.
2. DO NOT RSTERILIZE the ENDOTINE Triple™ device or its insertion tools.
3. Do not use if storage temperature exceeds 110°F (43°C) or the temperature indicator has turned black.
4. DO NOT cut or reconfigure the ENDOTINE Triple™ device.
5. Discard open, unused devices.
6. This device and its insertion tools are for SINGLE USE ONLY and MUST NEVER BE RE-USED.
7. Inspect all components for damage before use. Do not attempt to repair damaged components or sharpen the drill bit.
8. **Incomplete insertion of the device in the cranial bone may result in inadequate anchoring.**
9. Patients should be instructed to avoid contact sports or trauma to the head during healing.
10. Contour irregularities of the skin may be noticeable overlying the ENDOTINE Triple™ device, especially in non hair-bearing areas.

ADVERSE EFFECTS

Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissues can result in histological reactions.

Caution: Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

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Coapt Systems, Inc.