
Coapt™

ENDOTINE™ Forehead Device Instructions for Use

- English
Dansk
Deutsch
Español
Français
Italiano
Nederlands
Português
Svenska
Polski
Magyar
Ελληνικά
Русский
- CE 0344

ENGLISH

ENDOTINE Forehead™ Device Instructions for Use

DESCRIPTION
The ENDOTINE Forehead™ Device is a bioabsorbable fixation implant. This device is SINGLE USE ONLY and is supplied sterile. An instrument kit specific for this device is packaged and supplied separately.

INDICATIONS
The ENDOTINE Forehead™ Device is intended for use in subperiosteal browplasty surgery. The ENDOTINE Forehead™ Device is specifically indicated for use to fixate the sub-dermis to the cranial bone in browplasty.

CONTRAINDICATIONS
1. Patients appearing to have very thin bony thickness which might imply an inadequate cranial thickness for the ENDOTINE Forehead™ Device bone post which extends to a depth of 3.75 mm. The drilled hole extends to a depth of 3.95 mm +0.05, -0.0 (max 4.00 mm).
2. Situations where internal fixation is otherwise contraindicated (e.g. infection).
3. Thin, atrophic scalp.
4. Any known allergy or foreign-body sensitivities to plastic biomaterial, such as PLA/PGA.

LABELING, PACKAGING AND STERILIZATION
The ENDOTINE Forehead™ Device is sterilized by gamma irradiation. **Do NOT sterilize this device. Implant should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of punctures, tampering, water contamination or other damage.**

STORAGE INSTRUCTIONS
1. Store at room temperature (15 to 24 °C or 60 to 75 °F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the label.

INSTRUCTIONS FOR USE (for endoscopic or open techniques)
1. Prior to implant use, verify intraoperatively that the scalp thickness is appropriate for the ENDOTINE Forehead™ implants.
2. Make incisions large enough to allow insertion of the ENDOTINE Forehead™ Device without bending during placement; approximately 20 mm.
3. Perform the dissection in the usual fashion, making sure to achieve adequate release.
4. To minimize the risk of creating a full thickness cranial hole with the ENDOTINE Forehead™ drill bit, verify that the hole position is medial to the temporal fusion line and anterior to the coronal suture (Figure 1). Avoid the sagittal midline.
5. Plan the cranial hole placement such that the implant will ultimately lie under hair bearing tissue. Typically, for sagittal, paramedian incisions, the implant lies somewhat lateral and/or anterior to the incision. For coronal or transverse incisions, the implant lies anterior to the incision.
6. Prepare the cranial bone hole using a Coapt ENDOTINE Forehead™ drill bit (Figure 2) applied perpendicular to the bone. **After drilling, check that no residual loose bone material is present in the drill hole.**
7. Remove the ENDOTINE Forehead™ Device from the package foam using the Coapt Systems ENDOTINE Forehead™ Insertion Tool. Properly align the insertion tool with the ENDOTINE Forehead™ Device to avoid damage (Figure 3).
8. Utilizing the insertion tool, align the post of the device perpendicular to and within the center of the hole before pushing straight down (Figure 4).
9. With moderate pressure, insert the ENDOTINE Forehead™ Device into the prepared site until the device platform is resting on the cranium (Figure 5).
10. Release the ENDOTINE Forehead™ Device from the Insertion Tool. The platform of the ENDOTINE Forehead™ Device should rest flush against the cranium. **The device cannot be removed from the hole and reused.**
11. Perform the brow elevation and apply manual pressure to the scalp to ensure penetration of the tissue by the device (Figure 6).
12. If a drain is used, manipulate the subcutaneous portion so that it is maximal distance from the implant.
13. Use caution during skin closure to avoid elevating the scalp tissue off the implant.
14. A gentle pressure dressing may be appropriate upon completion of the case.

PRECAUTIONS
1. The implant ideally lies under hair bearing scalp at the conclusion of the procedure. Placement under forehead or alopecic skin may lead to unacceptable visibility of the device prior to absorption.

2. The patient is to be informed of possible adverse effects including palpability of the device prior to resorption, discomfort, and surgical risks. Tenderness and visible fullness over the device may be noted, especially in patients with thin scalp. These symptoms and/or wound problems may lead to a decision to surgically remove the device prior to full resorption.
3. Local cyst and/or sinus formation may occur and may be treated by aspiration or device removal.
4. High speed drilling (approximately 1000 rpm and greater) and repeat drilling may create holes that are wider than intended for a tight grip of the ENDOTINE Forehead™ Device. Single pass, low speed drilling reduces the occurrence of this issue.
5. Avoid any movement of the scalp, which may dislodge the tissue from the tines.
6. Support the head when turning or moving the patient.

WARNINGS
1. A drill bit is supplied for use in full size adult patients when used to create a hole in the cranium above (medial to) the temporal fusion line and anterior to the coronal suture. However, the surgeon should note any penetration beyond the inner table of the cranium and, if observed, should seek follow-up and neurological consultation as appropriate.
2. Do NOT RESTERILIZE the ENDOTINE Forehead™ Device.
3. Do not use if storage temperature exceeds 110 °F (43 °C) or the temperature indicator has turned black.
4. Do not cut or reconfigure the ENDOTINE Forehead™ Device.
5. Discard open, unused devices.
6. This device is for SINGLE USE ONLY and MUST NEVER BE RE-USED.
7. Visually inspect the insertion tool and drill bits for erosion, cracks or damage after each sterilization and prior to use. Do not attempt to repair or sharpen instruments.
8. Incomplete insertion of the device in the cranial bone may result in inadequate anchoring.
9. Patients should be instructed to avoid contact sports or trauma to the head during healing.
10. Contour irregularities of the skin may be noticeable overlying the ENDOTINE Forehead™ Device, especially in non hair-bearing areas.

ADVERSE EFFECTS
Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissues can result in histological reactions.

Caution: Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

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