Coapt

ENDOTINE Ribbon™





EMERGO EUROPE

Molenstraat 15
2513 BH, The Hague
The Netherlands

Coapt Systems, Inc.
1820 Embarcadero Road
Palo Alto, CA 94303

800.963.7670 Toll-free
650.461.7600 Main

0344

ENDOTINE Ribbon™

DESCRIPTION
The ENDOTINE Ribbon™ consists of a bioabsorbable implant supplied sterile for single use only.

INDICATIONS
The ENDOTINE Ribbon™ is indicated for use in elevation and fixation of tissues in the temporal region, midface, lower face/jowl to the deep temporal fascia and the neck to the mastoid fascia during cosmetic procedures (**Figure 1**).

CONTRAINDICATIONS

- Situations where internal fixation is otherwise contraindicated (e.g., infection).
- Patients appearing to have very thin soft tissues of the face in which the implant may be visible or palpable.
- Any known allergy or foreign-body sensitivities to plastic biomaterial.

PACKAGING, LABELING AND STERILIZATION
The ENDOTINE Ribbon™ is pre-sterilized by gamma irradiation and intended for single use only. Do **NOT** resterilize this device or the insertion tool. The device and its insertion tool should be accepted only if the factory packaging and labeling arrive intact. Do not use if packaging shows evidence of puncture(s), tampering, water contamination or other damage.

STORAGE INSTRUCTIONS
1. Store at room temperature (15 to 24°C or 60 to 75°F) in a dry place out of direct sunlight. Do not use beyond the expiration date listed on the outer label.

INSTRUCTIONS FOR USE

- Figure 2** identifies the fixation region (**2A**), anchoring leash (**2B**), and integrated release tab (**2C**).
- The implant is packaged within the insertion tool (**Figure 3**), a protective cover (**3A**) through which the release tab (**3B**) protrudes.
- The dissection (temporal, midface, lower face/jowl, neck) is performed using the surgeon's typical technique (see appropriate Dissection Plane listed in **Table A**).

Indication	Dissection Plane	Time Orientation
Temporal	Subcutaneous	Outward
Midface	Subperiosteal	Outward
Lower face/jowl	Subcutaneous, times engaging SMAS fascia	Inward
Neck	Superficial to platysmal muscle	Inward

Table A

- Insert the implant (tapered end first) through the incision into the desired position. If resistance is encountered, do not force the assembly into position, which may weaken the implant causing it to break. Remove the assembly and extend or enlarge your plane of dissection, with care to protect surrounding nerve structures, in order to provide adequate clearance for proper positioning of the device at the inferior recess of the dissected tissue.

ENGLISH

- Inform the patient of possible adverse effects including palpability of the device prior to resorption, discomfort, visibility of the device under the skin, and surgical risks inherent to the procedure (to include nerve injury, facial musculature dysfunction, infection, scarring, device erosion through the skin, asymmetry, loss of suspension/elevation, etc.). Some devices may require removal prior to their absorption due to discomfort, infection, reaction or other concerns. Some patients may desire temporary fillers to camouflage the borders of the implant during the healing phase.
- Transient local fluid accumulation and/or sinus formation may occur. Aspiration may yield traces of implant material and may result in healing of the sinuses without adverse effects.
- Use care around existing implanted devices.

WARNINGS

- DO NOT RSTERILIZE the ENDOTINE Ribbon™ or its insertion tool. This device and cover are for SINGLE USE ONLY and MUST NEVER BE RE-USED.
- Do not use if storage temperature exceeds 110°F (43°C) or the temperature indicator on the outer box has turned black.
- Discard open, unused devices.
- Exercise caution as the ENDOTINE Ribbon™ device contains several sharp tines that may catch or tear glows.
- Inspect all components for damage before use. Do not attempt to repair damaged components.

ADVERSE EFFECTS

Material sensitivity/allergic reactions in patients following surgery should be reported. Implantation of foreign materials in tissues can result in histological reactions.

Caution: Federal Law (USA) restricts this device to sale, distribution, or use by or on the order of a physician.

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